

Contact Day 1 Easy Wear

Technical Data

- Material: Methafilcon A
- Water content: 55 %
- Dk/t: 24
- Handling tint: light blue
- Central thickness (-3.00 D): 0.08 mm
- Front and Back surface: aspherical



ZEISS Contact Day 1 Easy Wear is a high-quality daily lens for spherical corrections.

Care recommendation

The lens does not need any care solution or product.


The characteristics

The Contact Day 1 Easy Wear is made for sporadic lens wearers and is characterised by the excellent comfort and visual acuity without any care expenditure. This is the reason why Contact Day 1 has established as a popular lens for regular lens wearer without using care products.

The geometry of the Contact Day 1 Easy Wear is characterized by an aspherical front and back surface and an extra slim edge design, whereby a very good centring is reached. Design and material properties make the handling of the Contact Day 1 very simple - also for untrained lens wearer.

The advantages of the Contact Day 1 Easy Wear

- High comfort
- Excellent visual acuity through aberration control
- Very good centring and flexibility
- Easy handling
- UV protection



Ocular Diagnostics

Zeiss Contact Lens Replacement Systems are distributed exclusively throughout Australia and New Zealand by

Ocular Diagnostics Australia Pty. Ltd.

ABN: 80 598 061 023

Mobile: 0423 203 245 All Hours
Phone/Facsimile: +61 (8) 83764193

Mail: P.O. Box 35
Glennelg SA 5045

Email: dallen@ocular diagnostics.com.au
Web: www.ocular diagnostics.com.au

Delivery program

Diameter	14.20 mm
Base curve	8.60 mm
Power (step)	-0.50 to -6.00 D (0.25 D)
	-6.50 to -9.00 D (0.50 D)
	+0.50 to +4.00 D (0.25 D)
	+4.50 to +6.50 D (0.50 D)
Packaging	30 lenses, 5 lenses