

Contact Day 30 toric



The Contact Day 30 toric is a back surface toric, dynamic stabilising monthly-replacement lens for astigmatic corrections.

Characteristics

The Contact Day 30 toric was designed for daily wear and offers the wearer with cylinder corrections the comfort and advantages of a modern replacement system.

The dynamic stabilisation prevents rotation of the lens on the cornea and holds the cylinder in the required axis position. This ensures visual acuity as well as sufficient mobility of the lens on the eye. The markings in 10° increments enable exact determination of the cylinder axis. As the lens is prism free, there are no vertical phorias with monocular use.

The thinness of this toric lens means that it is very oxygen permeable. The high-grade Ocufilcon D material with 55 % water content has good wetting properties and only a slight tendency to deposit formation. Design and material of the Contact Day 30 toric ensure outstanding wearing comfort and easy handling. Thanks to good centering, the Contact Day 30 toric contact lens is very easy to fit.

Advantages of the Contact Day 30 toric

- Excellent wearing comfort
- Reliable visual acuity and stabilisation
- Easy fitting
- Easy handling
- Light blue handling tint

Technical data

Material	Ocufilcon D with UV Protection
Water content	55 %
Dk	19×10^{-11}
Base curve	8.8 mm
Central thickness	0.12 mm (-3.00 D)
Geometry	spheric front surface, toric back surface

Delivery program

- Diameter: 14.40 mm
- Radius: 8.80 mm
- Sphere: plano to -6.00 D (increments 0.25)
- Cylinder: -0.75 D, -1.25 D, -1.75 D
- Axis: 10° to 180° (increments 10°)
- Packaging: blister pack, 6 per box

Care recommendation

- All soft contact lens care solutions and systems can be used.

